


Sequencing-Based DST at National PHL DST Reference Center for *Mycobacterium tuberculosis* when pure culture submitted

It is the responsibility of the submitting PHL to report results to the Ordering Clinician and TB Control Program

WGS DST TAT: 10-21 days

WGS Genotyping: The National PHL DST Reference Center for MTB will upload WGS data to the CDC National Surveillance System




Sequencing-Based DST at National PHL DST Reference Center for *Mycobacterium tuberculosis* when sediment, processed tissue, sterile body fluid, contaminated or non-viable culture submitted

It is the responsibility of the submitting PHL to report results to the Ordering Clinician and TB Control Program

tNGS TAT: 7-10 days
WGS DST TAT: 10-21 days

WGS Genotyping:
The National PHL DST Reference Center for MTB will upload WGS data to the CDC National Surveillance System



pDST panel that is set up depends on order. PZA pDST only set up when WGS is U for PZA. See Flow-Charts for pDST for details.


First-Line DST (FL-DST) at National PHL DST Reference Center for *Mycobacterium Tuberculosis*

Submit pure culture on LJ Media or MGIT Broth
-OR-
Specimen Reflexed from tNGS

Submitting PHL should keep 1 mL aliquot and submit remainder to MDL Isolates on Solid Media Acceptable

It is the responsibility of the submitting PHL to report results to the Ordering Clinician and TB Control Program

First Line Panel*
Rifampin (RIF), 1 µg/mL
Isoniazid (INH), 0.1 µg/mL[^]
Ethambutol (EMB), 5 µg/mL
Purity plate also setup to check for contamination or NTM
*pDST for PZA is performed only if "Uncertain" WGS results are obtained.
[^]INH tested at 0.4 µg/mL when INH resistance detected by any method




No Resistance Detected (including genotypic data)

Resistance Detected to at least one first-line drug

If RIF-R is detected, the sample will be referred to CDC for MDDR
PZA mono-resistance reflexed to SL-DST by request only

MDL Reports Results to Submitter

Smear or purity plate checked to rule out contamination




If no mutation was detected or amplification fails results will be reported.
No further action will be taken pending receipt of pure culture

DST for the 4-month regimen at National PHL DST Reference Center for *Mycobacterium Tuberculosis*

Submit pure culture on LJ Media or MGIT Broth
-OR-
Specimen Reflexed from tNGS

Submitting PHL should keep 1 mL aliquot and submit remainder to MDL Isolates on Solid Media Acceptable

It is the responsibility of the submitting PHL to report results to the Ordering Clinician and TB Control Program



MIGT Panel*
Rifampin (RIF), 1 µg/mL
Isoniazid (INH), 0.1 µg/mL[^]
Moxifloxacin (MOX), 0.25 µg/mL
Purity plate also setup to check for contamination or NTM
*pDST for PZA is set up only if "Uncertain" WGS results are obtained.
[^]INH tested at 0.4 µg/mL when INH resistance detected by any method

No Resistance Detected (molecular and first line testing reported as results become available)

Resistance Detected to RIF, INH or MOX

SL-DST will include EMB instead of MOX
If RIF-R is detected, the sample will be referred to CDC for MDDR
PZA mono-resistance reflexed to SL-DST by request only




MDL Reports Results to Submitter

Smear or purity plate checked to rule out contamination

No Contamination Detected
Phenotypic and Genotypic Results Concordant

No Contamination Detected
Phenotypic and Genotypic Results Discrepant -
OR- No genotypic data available


Contamination Detected



If no mutation was detected or amplification fails results will be reported.
No further action will be taken pending receipt of pure culture

Second-Line DST (SL-DST) at National PHL DST Reference Center for *Mycobacterium Tuberculosis*

It is the responsibility of the submitting PHL to report results to the Ordering Clinician and TB Control Program



If RIF-R is detected, the sample will be referred to CDC for MDDR. PZA mono-resistant isolates will not be automatically reflexed to SL-DST



Assessed for Concordance as Shown pg 1 and Below



Second Line Panel
 Amikacin (AMK), 1.5 µg/mL
 Capreomycin (CAP), 3 µg/mL
 Ethionamide (ETA), 5 µg/mL
 Kanamycin (KAN), 3.5 µg/mL
 Rifabutin (RFB), 0.5 µg/mL
 Moxifloxacin (MOX), 0.25 µg/mL

