

# Gram Stain Bench Guide

A reference guide to Gram stain preparation, interpretation and troubleshooting



This resource was produced by the Association of Public Health Laboratories (APHL) with special thanks to:

- Jaye Boman, MT(AMT), MLT(ASCP)
- Caryn Ivanof, M(ASCP)<sup>CM</sup>
- Laboratory Services for Northern Light Health—Eastern Maine Medical Center
- Mary DeMartino, MLS(ASCP)<sup>SM</sup>
- Patricia Ross, PhD, MLS(ASCP)<sup>CM</sup>
- Benjamin Hutton



**View a digital version of this resource**

**Introduction ..... 4**

Intended Use ..... 4

Gram Stain Reagents ..... 5

**Procedures..... 6**

Written Step-by-step Guide..... 6

Visual Step-by-step Guide..... 8

Reading Stained Slides..... 9

**Typical Staining ..... 10**

Gram-positive Organisms .....10

Gram-negative Organisms.....14

**Troubleshooting ..... 16**

References for Atypical Staining .....16

**Atypical Staining ..... 21**

Gram-positive Organisms ..... 21

Gram-negative Organisms..... 22

**Additional Resources..... 24**

## Intended Use

The Gram Stain Bench Resource Guide is designed to assist laboratory scientists in the accurate preparation and analysis of Gram-stained slides including the Gram staining procedure, typical staining, troubleshooting chart and atypical staining examples. This document is intended for laboratorians of all skill levels and is meant to supplement formal training. It is not intended to replace institutional protocols.

### Chart Key

**ANO<sub>2</sub>** .....Anaerobic environment

**BAP** .....Blood Agar Plate

**CHOC** .....Chocolate Agar Plate

**CO<sub>2</sub>** .....5-10% Carbon Dioxide

**CVA** .....*Campylobacter* Blood Agar

**MAC** .....MacConkey Agar Plate

**N<sub>2</sub>** .....85% Nitrogen environment

**O<sub>2</sub>** .....Oxygen/ambient air

**SDA** .....Sabouraud Dextrose Agar

**TSB** .....Tryptic Soy Broth

# Gram Stain Reagents

| Step                                 | Reagent  | Action  | Bacterial Cell Color After Step  |
|--------------------------------------|--|---|--|
| <b>Primary Stain</b>                 | Crystal violet   | Crystal violet penetrates the peptidoglycan in the bacterial cell wall.   | Purple   |
| <b>Mordant</b>                       | Gram's Iodine  | Iodine binds to the crystal violet, forming a complex that facilitates binding to cell wall.  | Purple   |
| <b>Decolorizer</b>                   | Either: <ul style="list-style-type: none"> <li>• 95% Ethanol</li> <li>• 95% Ethanol/ Acetone Mix*</li> </ul> | <ul style="list-style-type: none"> <li>• Bacteria with a thick peptidoglycan layer will resist decolorizing and retain the crystal violet-iodine complex.</li> <li>• Bacteria with a thin peptidoglycan layer will lose the crystal violet-iodine complex.</li> </ul> | <ul style="list-style-type: none"> <li>• Purple (Gram-positive)</li> <li>• No color (Gram-negative)</li> </ul> |
| <b>Counterstain/ Secondary Stain</b> | Safranin   | Safranin will enter the cell walls which have lost the crystal violet-iodine complex.   | <ul style="list-style-type: none"> <li>• Purple (Gram-positive)</li> <li>• Pink (Gram-negative)</li> </ul>     |

*\*Follow your SOP for the specific formula used in your laboratory.*

## Written Step-by-step Guide


See "Visual Step-by-step Guide" on page 8 for corresponding images to these instructions.

| Step                    | Action  |
|-------------------------|---|
| 1. Slide preparation    | <ul style="list-style-type: none"> <li>• If not using slides with imprinted circles, outline an area of the slide with a wax pencil.</li> <li>• <b>From an agar plate:</b> Add a drop of sterile saline or water inside the designated area on the slide. Using one end of a sterile wooden applicator stick or disposable needle touch one single, isolated colony and emulsify it into the water drop on the slide.</li> <li>• <b>From a broth:</b> Use a small disposable loop or needle to transfer a small amount of the liquid to the designated area on the slide.</li> <li>• Set the slide aside to dry fully.</li> </ul> |
| 2. Heat or Methanol Fix | <ul style="list-style-type: none"> <li>• <b>Heat Fix:</b> Once the slide is fully dry, place the slide onto a slide warmer set to 37°C. The heat will fix the cells to the slide so they don't wash off during staining. If the slide has not fully dried prior to heat fixing there is a risk of disrupting the cell walls or aerosolizing the cells, which is a potential biosafety risk.</li> <li>• <b>Methanol Fix:</b> Once slide is fully dry, flood slide with methanol for 60 seconds. Rinse methanol from slide using a wash bottle or running water while holding the slide at an angle.</li> </ul>                     |
| 3. Primary stain        | Flood the designated area of the slide with crystal violet and leave on for 60 seconds.   |
| 4. Rinse                | Using a wash bottle or running water and while holding the slide at an angle, rinse the slide to remove the crystal violet stain.   |


| Step                | Action  |
|---------------------|---|
| 5. Fix with mordant | Flood the slide with iodine and leave on for 60 seconds.  |
| 6. Rinse            | Using a wash bottle or running water and while holding the slide at an angle, rinse the slide to remove the iodine.   |
| 7. Decolorize       | <ul style="list-style-type: none"> <li>• Decolorize the slide by applying decolorizer and almost immediately rinsing with water.</li> <li>• Repeat this process until the purple color stops coming off while holding the slide at an angle.</li> </ul> <p><b>Note:</b> <i>If this is repeated beyond the time necessary to remove the crystal violet stain there is a risk of over-decolorizing which will result in Gram-positive cells appearing Gram-negative or Gram-variable. Proper timing will depend on decolorizer formulation.</i></p> |
| 8. Counterstain     | Flood the designated area of the slide with safranin and leave on for 30-60 seconds.  |
| 9. Rinse            | Using a wash bottle or running water and while holding the slide at an angle, rinse the slide to remove the counterstain.   |
| 10. Blot dry        | Gently blot excess liquid from slide and allow to air dry.  |

## Visual Step-by-step Guide

1. Slide preparation




2. Fix




3. Primary stain




4. Rinse




5. Fix with mordant




6. Rinse




7. Decolorize




8. Counterstain



9. Rinse







10. Blot dry








| Step                                | Action   |
|-------------------------------------|--|
| 1. Scan on low power                | Scan the slide at 10x to identify a well-stained area.   |
| 2. Scan on high power               | Scan the slide at 100x to be sure the proper area is still in focus and visible.   |
| 3. Add oil                          | <ul style="list-style-type: none"> <li>• Turn the nosepiece so no objective lens is touching the slide.</li> <li>• Without moving the slide or the stage, add a drop of immersion oil to the slide.</li> </ul> |
| 4. Read on high power oil immersion | Observe the previously identified area of the stained slide with a 100X high power objective lens.   |









**If the image goes out of focus or the well stained area is lost:** Wipe the oil from the lens and drop back to 10X to find the area of interest. The imprinted circle on the slide or the wax pencil mark can serve as a guide.




| Characteristics   |                | Cocci in Clusters   |   | Cocci in Pairs and Chains  |   |
|-------------------|----------------|---|---|--|---|
| Organism          |                | <i>Staphylococcus aureus</i>  |   | <i>Streptococcus pyogenes</i> (Gp A)   |   |
| Growth Conditions | Media State    | Solid Media   | Liquid Media  | Solid Media  | Liquid Media  |
|                   | Media          | BAP   | TSB   | BAP  | Blood Culture Bottle  |
|                   | Environment    | O <sub>2</sub>  | O <sub>2</sub>  | O <sub>2</sub>   | O <sub>2</sub>  |
|                   | Temperature    | 37°C  | 37°C  | 37°C   | 37°C  |
|                   | Time to Growth | 24 Hours  | 18 Hours  | 24 Hours   | 24 Hours  |
| Staining          |                |              |  |    |  |
| Notes             |                | Typically appear in clusters, similar to a cluster of grapes; detail in inset. Stains purple. |   | Typically appear as chains, similar to a string of beads or pearls; better chaining morphology in broth. Detail in inset. Stains purple. |   |

| Characteristics   |  | Cocci in Pairs and Chains  |   |                 |
|-------------------|--|--|---|-----------------|
| Organism          |  | <i>Streptococcus pneumoniae</i>  |   |                 |
| Growth Conditions | Media State  | Solid Media  | Solid Media   | Liquid Media    |
|                   | Media  | BAP  | BAP   | TSB             |
|                   | Environment  | CO <sub>2</sub>  | CO <sub>2</sub>   | CO <sub>2</sub> |
|                   | Temperature  | 36°C   | 36°C  | 36°C            |
|                   | Time to Growth   | 18 Hours   | 24 Hours  | 24 Hours        |
| Staining          |   |      |  |                 |
| Notes             | Typically appear in pairs, almost oval with the pointier ends together (lancet-shaped); detail in inset.<br>Stains purple. | <b>Capsule</b> —appears as a clearing around cells; detail in inset.<br>Stains purple. | More defined morphology from broth; detail in inset.<br>Stains purple.              |                 |

| Characteristics   |                | Rods with Spores  |  |
|-------------------|----------------|---|--|
| Organism          |                | <i>Bacillus subtilis</i>  | <i>Corynebacterium diphtheria</i>  |
| Growth Conditions | Media State    | Solid Media   | Solid Media  |
|                   | Media          | BAP   | BAP  |
|                   | Environment    | O <sub>2</sub>  | O <sub>2</sub>   |
|                   | Temperature    | 36°C  | 37°C   |
|                   | Time to Growth | 18 Hours  | 24 Hours   |
| Staining          |                |                                        |  |
| Notes             |                | Look for spores and their placement within the cell (terminal, subterminal or central); detail in inset. Stains purple. | Look for palisading/picket fence arrangement; detail in inset. Stains purple.      |

| Characteristics   |                | Yeast   |  |   |
|-------------------|----------------|---|--|---|
| Organism          |                | <i>Cryptococcus</i> spp.  | <i>Candida tropicalis</i>  |   |
| Growth Conditions | Media State    | Solid Media   | Solid Media  | Liquid Broth  |
|                   | Media          | SDA   | SDA  | TSB   |
|                   | Environment    | CO <sub>2</sub>   | O <sub>2</sub>   | O <sub>2</sub>  |
|                   | Temperature    | 37°C  | 37°C   | 7°C   |
|                   | Time to Growth | 48 Hours  | 72 Hours   | 24 Hours  |
| Staining          |                |  |  |  |
| Notes             |                | See budding detail in insets.<br>Stains purple.                                   |  |   |

| Characteristics   |   | Rods   |   |                             |
|-------------------|---|--|---|-----------------------------|
| Organism          |   | <i>Escherichia coli</i>  |   | <i>Campylobacter jejuni</i> |
| Growth Conditions | Media State   | Solid Media  | Liquid Media  | Solid Media                 |
|                   | Media   | BAP  | TSB   | CVA                         |
|                   | Environment   | O <sub>2</sub>   | O <sub>2</sub>  | N <sub>2</sub>              |
|                   | Temperature   | 37°C   | 37°C  | 42°C                        |
|                   | Time to Growth  | 24 Hours   | 24 Hours  | 48 Hours                    |
| Staining          |  |  |  |                             |
| Notes             | Stains pink.  |  | Displays gull-wing and "s" shapes.<br>Stains pink.                                  |                             |

| Characteristics   |                | Rods  |  | Diplococci  |
|-------------------|----------------|---|--|---|
| Organism          |                | <i>Haemophilus influenzae</i>   | <i>Fusobacterium</i> spp.  | <i>Neisseria gonorrhoeae</i>  |
| Growth Conditions | Media State    | Solid Media   | Solid Media  | Solid Media   |
|                   | Media          | CHOC  | BAP  | CHOC  |
|                   | Environment    | CO <sub>2</sub>   | ANO <sub>2</sub>   | O <sub>2</sub>  |
|                   | Temperature    | 37°C  | 37°C   | 37°C  |
|                   | Time to Growth | 24 Hours  | 48 Hours   | 24 Hours  |
| Staining          |                |  |  |    |
| Notes             |                | Displays pleomorphic morphology (coccobacilli).<br>Stains pink.                   | Cells have pointed ends.<br>Stains pink.   | Displays a “coffee or kidney bean” arrangement with the long, curved edges together in each cell of the pair.<br>Stains pink. |

This section outlines common issues that arise during Gram staining, along with possible causes and remediation techniques. The issue numbers are referenced in the Atypical Staining section.

## 1. Gram-positive appears pink/morphology inconsistent with Gram stain reaction

*e.g., Gram-negative cocci in chains*

| Problem                 | Possible Causes  | Remediation Options   |
|-------------------------|------------------|---|
| Technique               | Over decolorized | <ul style="list-style-type: none"> <li>• Verify decolorizing time</li> <li>• Use a timer</li> <li>• Consult SOP to confirm proper technique</li> <li>• Ensure acetone/alcohol ratio of decolorizer is appropriate for procedure times (the higher percentage of acetone, the shorter the decolorization time)</li> <li>• Ensure proper timing of exposure to iodine</li> <li>• Concurrently stain a QC slide</li> </ul> |
| Growth issues/<br>other | Culture is old   | <ul style="list-style-type: none"> <li>• Check date of subculture (encourage dating of media when inoculated)</li> <li>• Use fresh bacterial growth</li> </ul>  |

## 2. Gram-negative appears purple/morphology inconsistent with Gram stain reaction

| Problem                 | Possible Causes    | Remediation Options   |
|-------------------------|--------------------|---|
| Technique               | Smear is too thick | Prepare smear using a small drop of water and a loopful of organism   |
|                         | Under decolorized  | Verify decolorizing time  |
| Growth issues/<br>other | Culture is old     | <ul style="list-style-type: none"><li>• Check date of subculture (encourage dating of media when inoculated)</li><li>• Use fresh bacterial growth</li></ul> |

## 3. Inorganic material on slide/unable to see organisms or see pieces of organisms

| Problem  | Possible Causes       | Remediation Options  |
|--|-----------------------|--|
| External factors<br>(environment/<br>reagents) | Stain precipitate     | Filter stains, particularly crystal violet   |
|  | Slide is contaminated | <ul style="list-style-type: none"><li>• Keep box of slides covered when not in use</li><li>• Avoid potential splashing when making slides</li><li>• Ensure rinsing does not run over on other slides</li></ul> |

| Problem                         | Possible Causes   | Remediation Options   |
|---------------------------------|---|---|
| <b>External factors, cont'd</b> | Old staining reagents                                   | Use in-date reagents  |
|                                 | Improperly prepared/stored reagents                     | <ul style="list-style-type: none"> <li>• Confirm proper reagent preparation, if made in-house</li> <li>• Confirm reagents are stored properly</li> </ul>  |
|                                 | Immersion oil is contaminated                           | <ul style="list-style-type: none"> <li>• Confirm sterility of immersion oil/open new reagent</li> <li>• Clean objective</li> <li>• Review proper application of oil—do not touch slide</li> </ul> |
| <b>Technique</b>                | Staining slide before it is dry                         | Ensure slide is fully air-dried prior to fixation   |
|                                 | Excessive heat fixation                                 | Ensure slide is not held in heat too long (disrupts cell walls)   |
| <b>Growth issues/ other</b>     | The patient may be undergoing antimicrobial therapy     | If possible, confirm sample was collected prior to antibiotic treatment   |
|                                 | Organism not cultured in ideal growth environment/media | Confirm organism is growing in appropriate growth environment   |
|                                 | Organism too old/cell wall compromised                  | Subculture organism to confirm viability/stain fresh growth   |

#### 4. Unusual morphology/unable to determine morphology

| Problem                 | Possible Causes   | Remediation Options   |
|-------------------------|---|---|
| Technique               | Staining slide before it is dry                         | Ensure slide is fully air-dried prior to fixation                       |
|                         | Excessive heat fixation                                 | Ensure slide is not held in heat too long (disrupts cell walls)         |
| Growth issues/<br>other | The patient may be undergoing antimicrobial therapy     | If possible, confirm sample was collected prior to antibiotic treatment |
|                         | Organism not cultured in ideal growth environment/media | Confirm organism is growing in appropriate growth environment           |
|                         | Organism too old/cell wall compromised                  | Subculture organism to confirm viability/stain fresh growth             |




#### 5. Unable to see individual cells




| Problem   | Possible Causes   | Remediation Options   |
|-----------|---|---|
| Technique | <ul style="list-style-type: none"> <li>Smear too thick</li> <li>Smear too thin, prepared with too much water/not enough organism</li> </ul> | Prepare smear using a small drop of water and a loopful of organism |
|           | Focused on heavy part of smear  | Move microscope stage to thinner section of smear                   |




| Problem                 | Possible Causes   | Remediation Options   |
|-------------------------|---|---|
| Technique,<br>cont'd.   | Smear not fixed   | Fix organism to slide before staining                               |
|                         | Organism washed away  | Ensure proper rinsing/washing technique                             |
|                         | Smear area not clearly delineated                               | Delineate reading area with a wax pencil to help find correct plane |
|                         | Slide is upside down (organism facing down) on microscope stage | Flip slide on microscope stage                                      |
| Growth issues/<br>other | Organism too old/cell wall compromised                          | Subculture organism to confirm viability/stain fresh growth         |

## 6. Organism appears out of focus/unable to “zoom” in

| Problem                 | Possible Causes   | Remediation Options   |
|-------------------------|---|---|
| Technique               | Slide is upside down (organism facing down) on microscope stage | Flip slide on microscope stage                              |
|                         | No immersion oil used   | Add immersion oil   |
| Growth issues/<br>other | Organism too old/cell wall compromised                          | Subculture organism to confirm viability/stain fresh growth |

| Characteristics   |                             | Gram-positive Cocci   |   |   |
|-------------------|-----------------------------|---|---|---|
| Issue             | Problem                     | Over-decolorized  | Too Old   | Antibiotic-treated  |
|                   | Troubleshooting Reference # | 1   | 1, 2, 3, 5  | 3, 4  |
| Organism          |                             | <i>Staphylococcus aureus</i>  |   |   |
| Growth Conditions | Media                       | BAP   | BAP   | BAP   |
|                   | Environment                 | O <sub>2</sub>  | O <sub>2</sub>  | O <sub>2</sub>  |
|                   | Temperature                 | 37°C  | 37°C  | 37°C  |
|                   | Time to Growth              | 24 Hours  | 3 Days  | 24 Hours  |
| Staining          |                             |  |                                     |  |
| Notes             |                             | Cells can look Gram negative  | Cells can look inconsistently decolorized; they can appear stretched out or scant cells may be seen due to cell death | Displays bulbous/irregular shapes and scant number of cells                         |

| Characteristics   |                             | Gram-negative Rods  |   |   |
|-------------------|-----------------------------|---|---|---|
| Issue             | Problem                     | Over-decolorized  | Too Old   | Antibiotic-treated  |
|                   | Troubleshooting Reference # | 2   | 1, 2, 3, 5  | 3, 4  |
| Organism          |                             | <i>E. coli</i>  |   |   |
| Growth Conditions | Media                       | BAP   | BAP   | BAP   |
|                   | Environment                 | O <sub>2</sub>  | O <sub>2</sub>  | O <sub>2</sub>  |
|                   | Temperature                 | 37°C  | 37°C  | 37°C  |
|                   | Time to Growth              | 24 Hours  | 7 Days  | 24 Hours  |
| Staining          |                             |  |                                     |  |
| Notes             |                             | Cells can look Gram negative  | Cells can look inconsistently decolorized; they can appear stretched out or scant cells may be seen due to cell death | Displays bulbous/irregular shapes and scant number of cells                         |

| Characteristics   |                             | Gram-negative Rods  |  |   |
|-------------------|-----------------------------|---|--|---|
| Issue             | Problem                     | Smear Too Thick   | Stain Precipitate  | Smear Too Thick   |
|                   | Troubleshooting Reference # | 2, 5  | 3  | 2, 5  |
| Organism          |                             | <i>E. coli</i>  |  | <i>Haemophilus influenzae</i>   |
| Growth Conditions | Media                       | BAP   | BAP  | CHOC  |
|                   | Environment                 | O <sub>2</sub>  | O <sub>2</sub>   | CO <sub>2</sub>   |
|                   | Temperature                 | 37°C  | 37°C   | 36°C  |
|                   | Time to Growth              | 24 Hours  | 24 Hours   | 18 Hours  |
| Staining          |                             |  |  |  |
| Notes             |                             | Displays heavy section that looks Gram positive.                                  | Displays excess "blobs" of Crystal Violet  | Displays darker staining than usual, making it difficult to see morphology.         |

The following are additional helpful resources when determining best practices for performing a Gram stain:



## [Gram Stain Protocols](#)

American Society for Microbiology (ASM)



## [How to Perform a Gram Stain](#)

US Centers for Disease Control and Prevention (CDC)



## [Gram Stain Workshop for the Laboratory Generalist](#)

Nebraska Public Health Laboratory



## [The Gram Stain](#)

Virtual Interactive Bacteriology Laboratory, Michigan State University





## Association of Public Health Laboratories

The Association of Public Health Laboratories (APHL) works to strengthen laboratory systems serving the public's health in the US and globally. APHL's member laboratories protect the public's health by monitoring and detecting infectious and foodborne diseases, environmental contaminants, terrorist agents, genetic disorders in newborns and other health threats.

**7700 Wisconsin Avenue, Suite 1000 Bethesda, MD 20814 | 240.485.2745 | [www.aphl.org](http://www.aphl.org)**

© Copyright 2026, Association of Public Health Laboratories. All Rights Reserved.

This project was 100% funded with federal funds from a federal program of \$726,412. This publication was supported by Cooperative Agreement #NU600E000103 from the US Centers for Disease Control and Prevention (CDC). Its contents are solely the responsibility of the authors and do not necessarily represent the official views of CDC.